

January 2009

## **New Product Announcement Excyte Mini Automated ESR Analyzer – Code # EX-10310**

Inter Medico and Vital Diagnostics, Inc. are pleased to announce the latest addition to the Excyte ESR product line – the new Excyte Mini.

The Excyte Mini is an automated, 10 position, random access ESR analyzer, and it is perfect for low volume hospital laboratories. It has the same great features as the rest of the Excyte line including 10 minute pre-indication with results in 30 minutes, random access, automatic temperature correction, and a small footprint. Also available are the optional printer and barcode scanner.

The Excyte Mini uses the same 1.0 mL Excyte ESR tubes, in both vacuum and non vacuum configurations. This closed tube sampling and testing increases the safety of laboratory personnel by reducing exposure to blood.

The highly successful Excyte Placement Program is now available for the Excyte Mini system. In exchange for a three year commitment of Excyte ESR tubes and Accused ESR controls, customer will receive the Excyte Mini at no charge. System warranty and depot service are included free of charge for the duration of the contract term.

To obtain more information on the complete ESR product line, contact your sales representative at 1-800-387-9643 or at [info@intermedico.com](mailto:info@intermedico.com).

**excyte® mini**  
AUTOMATED ESR ANALYZER

